


Gentle Stim Control Foam

Neurostimulation Electrodes

- Clinically designed for all neurostimulation modalities like TENS and muscle stimulation
- Patented current controlling technology provides safe and comfortable electrical stimulation
- Optimized silver ink pattern printed onto conductive carbon film with soft white foam backing
- Hydrogel adhesive is optimized for maximum adhesion, durability and multiple applications to the skin
- Conductive pigtail wire is flexible and durable allowing for multiple uses
- Advanced design provides superior conformability to body and contours
- Available in multiple shapes and sizes
- Latex-free



Item #	Description	UOM
P850089	Foam, 1 1/4", Round	4/pk, 10 pk/bg
P850091	Foam, 2", Round	4/pk, 10 pk/bg
P850093	Foam, 2 3/4", Round	4/pk, 10 pk/bg
P850090	Foam, 2" x 2", Square	4/pk, 10 pk/bg
P850096	Foam, 2" x 4", Rectangle	4/pk, 10 pk/bg
P850098	Foam, 2" x 4", Oval	4/pk, 10 pk/bg

Distributed by:

Gentle Stim Iontophoresis

Electrode Treatment Kits

An innovative disposable iontophoretic drug delivery system, with a unique scientific blend of materials for optimum comfort during treatment.


- Each kit contains: one active drug delivery electrode, one buffered return electrode and two alcohol prep pads.

Active Drug Delivery Electrode

- Contains a reservoir pad to hold the deliverable solution.
- Reservoir pad fill area is conveniently accessible to allow for easy application of the solution prior to removal of the adhesive backing.
- The soft foam backing, with medical grade adhesive, is flexible and conforms to body contours.

Buffered Return Electrode

- The return electrode, with patented buffered gel, minimizes skin irritation.



Item #	Description	UOM
P850275	Small Square Electrode Kit <i>Includes: Small Square (2 3/4" x 2 3/4") Active Drug Delivery Electrode (1.5-2.0cc fill volume), 1 buffer return electrode, and 2 sterile, medium alcohol prep pads</i>	10 kits/bg
P850325	Medium Square Electrode Kit <i>Medium Square (3 1/4" x 3 1/4") Active Drug Delivery Electrode (2.5-3.0cc fill volume), 1 buffer return electrode, and 2 sterile, medium alcohol prep pads, 10 kt/bg</i>	10 kits/bg
P850350	Large Square Electrode Kit <i>Large Square (3 1/2" x 3 1/2") Active Drug Delivery Electrode (4.0-4.5cc fill volume), 1 buffer return electrode, and 2 sterile, medium alcohol prep pads, 10 kt/bg</i>	10 kits/bg
P850375	Butterfly Electrode Kit <i>Includes: Butterfly (3 3/8" x 3.38") Active Drug Delivery Electrode (2.0-2.5cc fill volume), 1 buffer return electrode, and 2 sterile, medium alcohol prep pads, 10 kt/bg</i>	10 kits/bg

Distributed by: