

Pro Advantage Lab Products

Pro
ADVANTAGE[®]

by **NDC**

Pregnancy Test Devices

- >99% Sensitivity; >99% Specificity
- 25 mIU/mL minimum hCG detection level
- Accurate results in 3 minutes for urine and 5 minutes for serum
- CLIA waived for urine



For Urine Samples,
Serum: Moderately Complex
CPT Code:
81025 Urine, 84703 Serum

Item #	Description
P080063	Urine hCG Pregnancy Cassette Device <i>Includes: 25 individually packaged devices and droppers. 25 tests per box</i>
P080021	Urine/Serum hCG Pregnancy Cassette Device <i>Includes: 25 individually packaged devices and droppers. 25 tests per box</i>
P080023	Urine hCG Pregnancy Test Strips <i>Includes: 20 Test Strips per canister, 5 canisters per box</i>

Strep A Test

- 97% Sensitivity; 95% Specificity
- Accurate results in 5 minutes
- CLIA waived and FDA 510(k) cleared to market
- Room temperature storage



CPT Code: 87880QW

Item #	Description
P080033	Strep A Test-Throat Swab Kit <i>Includes: 25 test strips, 25 swabs, 25 extraction tubes, 1-10mL reagent A, 1-10mL reagent B, 1-1mL positive control, 1-1mL negative control, package insert, procedure card. 25 tests per box</i>

Infectious Mononucleosis Test Device

- Qualitative membrane strip based immunoassay for detection of IM heterophile antibodies
- Simple to perform using whole blood (from venipuncture or fingerstick), serum or plasma
- Reliable results in 5 minutes
- CLIA waived for whole blood



For Whole Blood
Serum/Plasma:
Moderately Complex
CPT Code:
86308QW Whole Blood,
86308 Serum/Plasma

Item #	Description
P080016	Infectious Mononucleosis Test Device <i>Includes: 20 individually packaged devices, 20 disposable sample droppers, 20 disposable heparinized capillary tubes, 1 dispensing bulb, 1 sample buffer, 1 procedure card, 1 package insert. 20 tests per box</i>

Fecal Occult Blood Test Devices

- On board positive and negative monitors
- Includes enhanced color sensitivity and readability
- CLIA waived



CPT Code: 82270

Item #	Description
P080018	Fecal Occult Blood Lab Pack – Guaiac Slides <i>Includes: 100 Single Slides Tests, 2x10mL Developers, Applicators. 100 per box</i>
P080028	Fecal Occult Developer Solution – for Guaiac Slides <i>Includes: 8x10mL Developers. 8 per box</i>



H. pylori Test Device

- 88% sensitivity; 89% specificity
- Accurate results in 10 minutes
- CLIA waived for whole blood

Item #	Description
P080095	H. pylori Test Device – Whole Blood/Serum/Plasma <i>Includes: 20 individually packaged devices, 20 disposable sample droppers, 20 disposable heparinized capillary tubes, 1 dispensing bulb, 1 sample buffer, 1 procedure card, 1 package insert. 20 tests per box</i>



For Whole Blood, Serum/Plasma: Moderately Complex

CPT Code:
86318QW Whole Blood,
86677 Serum/Plasma

Tourniquets

- Latex-Free
- Tear resistant
- Textured for non-slip performance
- Individually rolled and banded
- Easily tied and untied
- Non-sterile



Item #	Description	Size	Color	Bag
P929900	Latex-Free Tourniquet	1" x 18"	White	100

Alcohol Prep Pads

- Meets the highest bacteriological standards
- Contains 70% Isopropyl Alcohol

Item #	Description	Box
P902045	Medium Alcohol Prep Pad, 1³/₁₆" x 2⁵/₈"	200 pk
P902050	Large Alcohol Prep Pad, 1³/₄" x 3¹/₂"	200 pk



Microscope Slides

- Thin, flat and resistant to corrosion or fogging
- Precision ground edges on all four sides
- Slides measure 25mm x 75mm x 1mm

Item #	Description	Box
P460025	Plain Microscope Slides	72
P460125	Frosted Microscope Slides	72
P460122	Cover Glass (0.13 to 0.17mm thick)	1 ounce
P460222	Cover Glass (0.17 to 0.25mm thick)	1 ounce



Specimen Containers

- Latex Free
- Non-Toxic

Item #	Description	Box
Sterile		
P250400	4 oz. Specimen Container with screw on lid & label, packaged individually in Poly Pouch	100
P250410	4 oz. Specimen Container with screw on lid & Tamper Evident Label	100
Non-Sterile		
P250415	4 oz. Specimen Container with screw on lid	500



Cotton Tipped Applicators

- Available either non-sterile or sterile wrapped

Item #	Description	Box
75900	Applicator Stick, 6" x 1/2"	1000
76200	Cotton Tipped Applicator, 6" x 1/2", Wooden Shaft, Non-Sterile	10 bags of 100
76700	Cotton Tipped Applicator, 6" x 1/2", Wooden Shaft, Sterile	100 packs of 1
77600	Rayon Tipped Applicator, 8", Plastic Shaft, Non-Sterile	50



Professional Towels

- Poly backing prevents moisture penetration

Item #	Description	Size	Color	Case
P752023	2 Ply Tissue + Poly	13" x 18"	White	500



Cotton Balls

- Made from absorbent 100% bleached cotton
- Soft, compact and consistent in size and weight
- Packaged in bag with resealable zip closure for easy storage

Item #	Description	Size	Color	Bag
P159025	Cotton Balls	Medium	White	2000
P159045	Cotton Balls	Large	White	1000



Urine Reagent Strips

- **10 parameter Urine Reagent Strip**

For the semi-quantitative and qualitative detection of Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose in urine.

- **7 parameter Urine Reagent Strips**

For the semi-quantitative and qualitative detection of Leukocytes, Nitrite, Protein, pH, Blood, Ketone, Glucose in urine.

- **4 parameter Urine Reagent Strips**

For the semi-quantitative and qualitative detection of Leukocytes, Nitrite, Protein, Glucose in urine.

- **2 parameter Urine Reagent Strips**

For the semi-quantitative and qualitative detection of Glucose and Protein in urine.

Item #	Description	Bottle
P080010	10 parameter Urine Reagent Strips	100
P080007	7 parameter Urine Reagent Strips	100
P080014	4 parameter Urine Reagent Strips	100
P080012	2 parameter Urine Reagent Strips	100



For Urinalysis, Visual
CPT Code: 81002
For Urinalysis, Automated
CPT Code: 81003QW

Urine Analyzer

Portable urine chemistry analyzer exclusively for reading Pro Advantage Urine Reagent Strips.

- Process 120 tests/hour
- Large 2,000 patient memory storage
- EMR/LIS Connectivity
- Eliminates color interpretation variability
- Minimizes transcription errors with an easy-to-read printout
- Patient ID Capable
- LCD screen
- User friendly, portable and convenient

Item #	Description
P080000	Urine Analyzer
P080001	Barcode Reader
P080002	Liquid Controls <i>Contains: 1-5 mL Bottle Positive Control, 1-5 mL Bottle Negative Control</i>
P080004	Thermal Printer Paper, 4 roll/bx
P080005	Data Transfer Cable



For Urinalysis, Visual
CPT Code: 81002
For Urinalysis, Automated
CPT Code: 81003QW

Distributed By: