

# Pro Advantage Lab Products

**Pro**  
**ADVANTAGE**<sup>®</sup>

by **NDC**

## Pregnancy Test Devices

- >99% Sensitivity; >99% Specificity
- 25 mIU/mL minimum hCG detection level
- Accurate results in 3 minutes for urine and 5 minutes for serum
- CLIA waived for urine



Item #	Description
--------	-------------

<b>P080063</b>	<b>Urine hCG Pregnancy Cassette Device</b>
----------------	--

*Includes: 25 individually packaged devices and droppers. 25 tests per box*

<b>P080021</b>	<b>Urine/Serum hCG Pregnancy Cassette Device</b>
----------------	--

*Includes: 25 individually packaged devices and droppers. 25 tests per box*

<b>P080023</b>	<b>Urine hCG Pregnancy Test Strips</b>
----------------	--

*Includes: 20 Test Strips per canister, 5 canisters per box*



**For Urine Samples,  
Serum: Moderately Complex**  
CPT Code:  
81025 Urine, 84703 Serum

## Strep A Test

- 97% Sensitivity; 95% Specificity
- Accurate results in 5 minutes
- CLIA waived and FDA 510(k) cleared to market
- Room temperature storage



Item #	Description
--------	-------------

<b>P080033</b>	<b>Strep A Test-Throat Swab Kit</b>
----------------	-------------------------------------

*Includes: 25 test strips, 25 swabs, 25 extraction tubes, 1-10mL reagent A, 1-10mL reagent B, 1-1mL positive control, 1-1mL negative control, package insert, procedure card. 25 tests per box*



CPT Code: 87880QW

## Infectious Mononucleosis Test Device

- Qualitative membrane strip based immunoassay for detection of IM heterophile antibodies
- Simple to perform using whole blood (from venipuncture or fingerstick), serum or plasma
- Reliable results in 5 minutes
- CLIA waived for whole blood



Item #	Description
--------	-------------

<b>P080016</b>	<b>Infectious Mononucleosis Test Device</b>
----------------	---

*Includes: 20 individually packaged devices, 20 disposable sample droppers, 20 disposable heparinized capillary tubes, 1 dispensing bulb, 1 sample buffer, 1 procedure card, 1 package insert. 20 tests per box*



**For Whole Blood  
Serum/Plasma:  
Moderately Complex**  
CPT Code:  
86308QW Whole Blood,  
86308 Serum/Plasma

## Fecal Occult Blood Test Devices

- On board positive and negative monitors
- Includes enhanced color sensitivity and readability
- CLIA waived



CPT Code: 82270

Item #	Description
<b>P080018</b>	<b>Fecal Occult Blood Lab Pack – Guaiac Slides</b> <i>Includes: 100 Single Slides Tests, 2x10mL Developers, Applicators. 100 per box</i>
<b>P080028</b>	<b>Fecal Occult Developer Solution – for Guaiac Slides</b> <i>Includes: 8x10mL Developers. 8 per box</i>



## H. pylori Test Device

- 88% sensitivity; 89% specificity
- Accurate results in 10 minutes
- CLIA waived for whole blood

Item #	Description
<b>P080095</b>	<b>H. pylori Test Device – Whole Blood/Serum/Plasma</b> <i>Includes: 20 individually packaged devices, 20 disposable sample droppers, 20 disposable heparinized capillary tubes, 1 dispensing bulb, 1 sample buffer, 1 procedure card, 1 package insert. 20 tests per box</i>



**For Whole Blood, Serum/Plasma: Moderately Complex**

CPT Code:  
86318QW Whole Blood,  
86677 Serum/Plasma

## Tourniquets

- Latex-Free
- Tear resistant
- Textured for non-slip performance
- Individually rolled and banded
- Easily tied and untied
- Non-sterile



Item #	Description	Size	Color	Bag
<b>P929900</b>	<b>Latex-Free Tourniquet</b>	<b>1" x 18"</b>	<b>White</b>	<b>100</b>

## Alcohol Prep Pads

- Meets the highest bacteriological standards
- Contains 70% Isopropyl Alcohol

Item #	Description	Box
Sterile		
<b>P902045</b>	<b>Medium Alcohol Prep Pad, 1<sup>3</sup>/<sub>16</sub>" x 2<sup>5</sup>/<sub>8</sub>"</b>	<b>200 pk</b>
<b>P902050</b>	<b>Large Alcohol Prep Pad, 1<sup>3</sup>/<sub>4</sub>" x 3<sup>1</sup>/<sub>2</sub>"</b>	<b>200 pk</b>



## Microscope Slides

- Thin, flat and resistant to corrosion or fogging
- Precision ground edges on all four sides
- Slides measure 25mm x 75mm x 1mm

Item #	Description	Box
P460025	Plain Microscope Slides	72
P460125	Frosted Microscope Slides	72
P460122	Cover Glass (0.13 to 0.17mm thick)	1 ounce
P460222	Cover Glass (0.17 to 0.25mm thick)	1 ounce



## Specimen Containers

- Latex Free
- Non-Toxic

Item #	Description	Box
Sterile		
P250400	4 oz. Specimen Container with screw on lid & label, packaged individually in Poly Pouch	100
P250410	4 oz. Specimen Container with screw on lid & Tamper Evident Label	100
Non-Sterile		
P250415	4 oz. Specimen Container with screw on lid	500



## Cotton Tipped Applicators

- Available either non-sterile or sterile wrapped

Item #	Description	Box
75900	Applicator Stick, 6" x 1/2"	1000
76200	Cotton Tipped Applicator, 6" x 1/2", Wooden Shaft, Non-Sterile	10 bags of 100
76700	Cotton Tipped Applicator, 6" x 1/2", Wooden Shaft, Sterile	100 packs of 1
77600	Rayon Tipped Applicator, 8", Plastic Shaft, Non-Sterile	50



## Professional Towels

- Poly backing prevents moisture penetration

Item #	Description	Size	Color	Case
P752023	2 Ply Tissue + Poly	13" x 18"	White	500



## Cotton Balls

- Made from absorbent 100% bleached cotton
- Soft, compact and consistent in size and weight
- Packaged in bag with resealable zip closure for easy storage

Item #	Description	Size	Color	Bag
P159025	Cotton Balls	Medium	White	2000
P159045	Cotton Balls	Large	White	1000



## Urine Reagent Strips

- **10 parameter Urine Reagent Strip**

For the semi-quantitative and qualitative detection of Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose in urine.

- **7 parameter Urine Reagent Strips**

For the semi-quantitative and qualitative detection of Leukocytes, Nitrite, Protein, pH, Blood, Ketone, Glucose in urine.

- **4 parameter Urine Reagent Strips**

For the semi-quantitative and qualitative detection of Leukocytes, Nitrite, Protein, Glucose in urine.

- **2 parameter Urine Reagent Strips**

For the semi-quantitative and qualitative detection of Glucose and Protein in urine.

Item #	Description	Bottle
P080010	<b>10 parameter Urine Reagent Strips</b>	<b>100</b>
P080007	<b>7 parameter Urine Reagent Strips</b>	<b>100</b>
P080014	<b>4 parameter Urine Reagent Strips</b>	<b>100</b>
P080012	<b>2 parameter Urine Reagent Strips</b>	<b>100</b>



For Urinalysis, Visual  
CPT Code: 81002  
For Urinalysis, Automated  
CPT Code: 81003QW

## Urine Analyzer

Portable urine chemistry analyzer exclusively for reading Pro Advantage Urine Reagent Strips.

- Process 120 tests/hour
- Large 2,000 patient memory storage
- EMR/LIS Connectivity
- Eliminates color interpretation variability
- Minimizes transcription errors with an easy-to-read printout
- Patient ID Capable
- LCD screen
- User friendly, portable and convenient

Item #	Description
P080000	<b>Urine Analyzer</b>
P080001	<b>Barcode Reader</b>
P080002	<b>Liquid Controls</b> <i>Contains: 1-5 mL Bottle Positive Control, 1-5 mL Bottle Negative Control</i>
P080004	<b>Thermal Printer Paper, 4 roll/bx</b>
P080005	<b>Data Transfer Cable</b>



For Urinalysis, Visual  
CPT Code: 81002  
For Urinalysis, Automated  
CPT Code: 81003QW

Distributed By: